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Blood Glucose Monitoring System

User Manual

For single-patient use / home use only.
For use with ReliOn Premier Test Strips only.
Test strips sold separately.
Lancets and lancing device sold separately.

Welcome to the ReliOn™ Premier BLU Blood Glucose Monitoring System

Thank you for choosing the ReliOn™ Premier BLU Blood Glucose Monitoring System. The system provides you with rapid and convenient blood glucose *in vitro* (i.e., outside the body) diagnostic monitoring. You can obtain accurate results in just five seconds with a small (0.5 µL) blood sample.

Test results can be sent to a smartphone through Bluetooth® wireless connection.

Bluetooth® wireless technology

Bluetooth® wireless technology is used by some smartphones and many other devices. ReliOn™ Premier BLU Blood Glucose Meter uses Bluetooth® wireless technology to pair* and send your blood glucose test results to a smartphone.

The ReliOn Premier BLU Blood Glucose Meter is designed to work with a mobile app. When using the ReliOn™ Premier BLU Blood Glucose Monitoring System, we suggest you pair your ReliOn™ Premier BLU Blood Glucose Meter with your smartphone and track your blood glucose test results.

- No part of this document may be reproduced in any form or by any means without the prior written consent of the manufacturer.
- The information in this manual is correct at the time of printing. However, the manufacturer reserves the right to make any necessary changes at any time without notice as our policy is one of continuous improvement.

* The process of creating a connection between two Bluetooth® devices. An auto-generated passkey has to be exchanged between the two devices. Once the devices are paired (connected), they will automatically communicate with each other when the Bluetooth® feature is activated.

DISTRIBUTED BY:

Walmart Inc.,
Bentonville, AR 72716
www.ReliOnBGM.com

© All Rights Reserved.
P/N 5409-06 Rev. 08/24

PGUAA-0000368 REV0 2024-08

Customer Service

1-833-256-3561
Monday-Friday 7 a.m. to 7 p.m.
Central time. Weekends and Holidays
8 a.m. to 5 p.m. Central time.
In case of emergency, please contact
your healthcare professional or
emergency medical response.

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WARNING

Turn off the Bluetooth® feature in areas where the use of wireless devices is restricted, such as hospitals, airplanes, and some healthcare professional offices.

Trademarks

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth® SIG, Inc. and any use of such marks by the manufacturer is under license. All other trademarks and trade names are those of their respective owners.

FCC Radiation Exposure Statement

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment.

This equipment should be installed and operated with minimum distance of eight inches (20 cm) between the radiator and your body.

Any changes or modifications not expressly approved by the manufacturer could void the user's authority to operate this equipment.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Page intentionally left blank

Table of Contents

Information

Important Information: Read This First	6
Important Safety Information	7
Limitations of ReliOn™ Premier BLU Blood Glucose Monitoring System	8
Specifications	10
ReliOn™ Premier BLU Blood Glucose Monitoring System	12
Inserting or Replacing the Batteries	13
ReliOn™ Premier Blood Glucose Test Strip	14
ReliOn™ Premier BLU Blood Glucose Meter	17
ReliOn™ Premier BLU Blood Glucose Meter Display	18
Understanding the Color Screen	19

Preparation

Setting Up the ReliOn™ Premier BLU System	20
Setting Up Bluetooth®	21
Adjusting the Date and Time	23
Turning the Sound On/OFF	25
Turning on the Test Strip Expiration Date Indicator	26
Setting the Target Range	27
Target Blood Glucose Ranges	29
Setting the Test Strip Expiration Date Indicator	30
Checking the System	31
Control Solution Testing	32
Comparing the Control Solution Test Results	35

Testing

Using the Lancing Device	36
Preparing the Lancing Device	37
Preparing the Meter and Test Strip	38
Applying Blood Sample	39
Discarding Used Lancets	43
'Hi' and 'Lo' Messages	44
Transferring Test Results Using the Data Cable	45

Additional Functions

Meter Memory	46
Viewing Past Test Averages	46
Viewing Test Results	49
Setting the Alarm Function	50
Setting the Post-meal Alarm (PP2 alarm)	50
Setting the Time Alarms (Alarm 1–3)	52
Caring for the ReliOn™ Premier BLU System	54
Cleaning and Disinfection Procedure	56

Maintenance

Understanding Error Messages	60
General Troubleshooting	62
Performance Characteristics	63
EMC Table	65
Warranty Information	73

Important Information: Read This First

To receive safe and optimum system benefits, please read the entire manual contents before using the system. Please note the following information.

- **Do not** use the system for diagnosis or screening of diabetes.
- Use only fresh capillary whole blood samples for testing.

Intended use

- The ReliOn™ Premier BLU Blood Glucose Monitoring System is comprised of the ReliOn™ Premier BLU Blood Glucose Meter and the ReliOn™ Premier Blood Glucose Test Strips.
- The ReliOn™ Premier BLU Blood Glucose Monitoring System is intended for the quantitative measurement of glucose in fresh capillary whole blood samples drawn from the fingertips.
- The ReliOn™ Premier BLU Blood Glucose Monitoring System is intended for self-testing outside the body (*in vitro diagnostic use*) by people with diabetes at home as an aid to monitor the effectiveness of diabetes control.
- The system is intended to be used by a single person and should not be shared.
- The system is not intended for use on neonates and is not for the diagnosis or screening of diabetes.

Meaning of Symbols Used:

 Caution

 Biological risks

Important Safety Information

- The ReliOn™ Premier Control Solutions are for use with the ReliOn™ Premier BLU Blood Glucose Meter and ReliOn™ Premier Test Strips to check that the meter and test strips are working together properly, and that the test is performing correctly.
- Please use this device only for the intended use described in this user manual.
- Please follow the suggested cleaning and disinfection procedures described in this user manual.
- ReliOn™ Premier Blood Glucose Test Strips are intended for single use only. They should be disposed of in an appropriate container immediately after use.
- The ReliOn™ Premier BLU Blood Glucose Meter is designed to minimize code related errors in monitoring by using the no-coding function.
- Glucose in blood samples react with the chemical in the test strip to produce a small electrical current. The ReliOn™ Premier BLU Blood Glucose Meter detects this electrical current and measures the amount of glucose in the blood sample.
- The target range is set by default at 70–180 mg/dL. Consult your healthcare professional to activate the Target Range Indicator and change the target range for your needs.
- Use of accessories, transducers and cables other than those specified or provided by the authorized distributor could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Limitations of ReliOn™ Premier BLU Blood Glucose Monitoring System

- An abnormally high or low red blood cell count (hematocrit level over 65 % or below 15 %) may produce inaccurate results.
- Inaccurate results may occur in severely hypotensive individuals or patients in shock.
- Inaccurate low results may occur for individuals experiencing a hypoglycemic hyperosmolar state, with or without ketosis.
- Severe dehydration (excessive water loss) may cause false low results. If you believe you are suffering from severe dehydration, consult your healthcare professional immediately.
- Altitudes higher than 10,000 ft. (3,048 m) above sea level may have an effect on the performance of the test strip.
- This system is for single-patient use only and should not be shared.
- Please note that individuals with color vision deficiency (CVD) may experience difficulty with the Target Range Indicator of this product. This function can be turned on and off as needed, but users with CVD should be aware that it may not be applicable due to potential color differentiation issues. If you have any questions or concerns, please consult with a healthcare professional or contact our customer service team for assistance.
- **Not** for neonatal use.
- **Not for use** on critically ill patients.
- **Not for use** in anyone undergoing oxygen therapy.
- **For *in vitro* diagnostic use only.**

Limitations of ReliOn™ Premier BLU Blood Glucose Monitoring System

- **Do not use** for the diagnosis of or screening of diabetes mellitus.
- The ReliOn™ Premier BLU Blood Glucose Meter should only be used with ReliOn™ Premier Blood Glucose Test Strips.
- For questions or concerns, contact Customer Service at 1-833-256-3561 seven days a week. Standard hours are Monday-Friday 7 a.m. to 7 p.m. Central time. Saturday-Sunday (and holiday) hours are 8 a.m. to 5 p.m. Central time. At all other times or in case of emergency, please contact your healthcare professional or emergency medical response.
- Portable RF Communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 in) to any part of the ReliOn™ Premier BLU meter.
- Strong electromagnetic fields may interfere with the proper operation of the meter. Do not use the meter close to sources of strong electromagnetic radiation.
- To avoid electrostatic discharge, do not use the meter in a very dry environment, especially one in which synthetic materials are present.

This device is not intended for use in healthcare or assisted-use setting such as hospitals, physician offices or long-term care facilities because it has not been cleared by the FDA for use in these settings, including for routine assisted testing or as part of glycemic control procedures.

Use of this device on multiple patients may lead to the transmission of Human Immunodeficiency Virus (HIV), Hepatitis C Virus (HCV), Hepatitis B Virus (HBV), or other bloodborne pathogens.

Specifications

Product specifications	
Measurement range	20–600 mg/dL
Sample size	Minimum 0.5 µL
Test time	Five seconds
Sample type	Fresh capillary whole blood
Calibration	Plasma-equivalent
Assay method	Electrochemical
Battery life	1,000 tests
Power	Two 3.0 V lithium batteries (disposable, type CR2032)
Memory	1,000 test results
Size	4.06 x 2.23 x 0.66 inches (103 x 54 x 16.7 mm)
Weight	2.53 oz (71.7 g) with batteries
Bluetooth® technology	<ul style="list-style-type: none">Frequency range: 2.4–2.4835 GHzOperating range distance: maximum 32 feet (10 m) unobstructedOperating channels: 40 channelsSecurity encryption: 128-bit AES (Advanced encryption standard)

Operating ranges

Temperature	42.8–111.2 °F (6–44 °C)
Relative humidity	10–90 %
Hematocrit	15–65 %

Specifications

Storage/Transport conditions		
Temperature	Glucose Meter (with batteries)	32–122 °F (0–50 °C)
	Test strip	34–86 °F (1–30 °C)
	Control solution	46–86 °F (8–30 °C)
Relative humidity	Test strip	20–80 %

ReliOn™ Premier BLU Blood Glucose Monitoring System

ReliOn™ Premier BLU Blood Glucose Monitoring System includes the following items:

- ReliOn™ Premier BLU Blood Glucose Meter
- User Manual
- Quick Reference Guide
- Bluetooth® Pairing Quick Guide
- Logbook
- Carrying Case
- Batteries (2)

Check all the components after opening the ReliOn™ Premier BLU Blood Glucose Monitoring System package.

ReliOn™ Premier Test Strips and Control Solutions are necessary but not included.

Lancing device and lancets are needed for testing but not included.

ReliOn™ Premier Control Solutions and data cable for the diabetes management system can be ordered separately by contacting Customer Service at 1-833-256-3561. At all other times, contact your healthcare professional.

The lancets and lancing device can be purchased through online retailers (e.g. Walmart.com) or by telephone by contacting Customer Service at 1-833-256-3561.

CAUTION

- The meter and lancing device are for single-patient use.
- **DO NOT** share them with anyone including other family members.
- **DO NOT** use on multiple patients.

Inserting or Replacing the Batteries

- The ReliOn™ Premier BLU Blood Glucose Meter uses two 3.0 V lithium batteries.
- **Do not** mix old and new batteries.
- **Do not** mix alkaline, standard (carbon-zinc) or rechargeable (ni-cad, ni-mh, etc.) batteries.
- Before using the meter, check the battery compartment and insert batteries if empty.
- When the  symbol appears on the display while the meter is in use, the batteries should be replaced as soon as possible.
- The test results may not be saved if the batteries run out.

Step 1

Make sure the meter is turned off. Push the cover in the direction of the arrow to open the battery compartment.

Step 2

Remove the used batteries and insert two new batteries with the '+' side facing up. Make sure the batteries are inserted firmly.

Step 3

Place the cover on the battery compartment. Close the battery cover firmly. A click may be heard.



NOTE

- Removing the meter batteries will not affect your stored results. However, meter setting may need to be reset, see page 18.

ReliOn™ Premier Blood Glucose Test Strip

The ReliOn™ Premier BLU Blood Glucose Monitoring System measures blood glucose quickly and accurately. It automatically absorbs the small blood sample applied to the narrow edge of the test strip.

Contact bars

Gently push the test strip, with its contact bars facing up, into the test strip port of meter.

Confirmation window

Check here to see whether sufficient blood sample has been applied.

Edge to apply blood sample

Apply blood sample here for testing.

Warning

- ReliOn™ Premier Blood Glucose Test Strips should only be used with fresh capillary whole blood samples.
- **Do not** reuse test strips.
- **Do not** use test strips past the expiration date.
- When test strips are used and stored according to correct storage and handling methods, both new, unopened vials and vials that have been opened can be used until the expiration date printed on the test strip box and vial label.
- Store test strips in a cool and dry place at a temperature between 34–86 °F (1–30 °C) and 20–80 % relative humidity.
- Keep test strips away from direct sunlight or heat, and **do not** freeze.
- Store test strips only in their original vial.
- Close the vial tightly after taking out a test strip for testing and use the test strip immediately.
- Avoid getting any liquid or moisture in the test strip vial. This can affect the test strips and cause inaccurate test results.
- **Do not** apply samples other than capillary whole blood or control solution to the test strip.
- Handle test strips only with clean and dry hands.
- **Do not** bend, cut or alter test strips in any way.
- For detailed storage and usage information, refer to the ReliOn™ Premier Blood Glucose Test Strip package insert.

CAUTION

- Keep the meter and testing supplies away from young children.
- Drying agents in the vial cap may be harmful if inhaled or swallowed, and may cause skin or eye irritation.

CAUTION

- **Do not** allow any foreign substances or liquid substances, such as dirt, blood, or water, enter into the meter. The meter may be damaged or may malfunction. Follow the warning information provided below to prevent possible damage to the meter.
- **Do not** apply the blood or control solution samples directly to the test strip port.
- **Do not** apply the blood or control solution samples to the test strip while holding the meter in a way that the tip of the test strip faces upwards. The blood or control solution samples may run down the surface of the test strip and flow into the test strip port.
- **Do not** store your meter in unsanitary or contaminated sites.
- Make sure to follow the Cleaning and Disinfection Procedures found in the Caring for the ReliOn™ Premier BLU System section of this user manual.

NOTE

You can get more safety information at CDC Clinical Reminder (www.cdc.gov/injectionsafety/Fingerstick-DevicesBGM.html).

ReliOn™ Premier BLU Blood Glucose Meter

Data Port

Used to transfer data from the meter to a computer with a data cable.

<, > Button

Turns the meter on, selects or changes information.

Display

Shows results and messages.

○ Button

Turns the meter on/off, confirms menu selections, and changes information.

Test Strip Port

Insert test strip here.

Test Strip Ejector

Slide down to discard used test strips.

NOTE

- The data cable for the data management software can be ordered by contacting Customer Service at 1-833-256-3561.
- The unit of measurement is fixed and cannot be changed by the user.

ReliOn™ Premier BLU Blood Glucose Meter Display

- ⑥ **Test results:** test results display panel
- ⑦ **alarm:** appears when the time alarm has been set
- ⑧ **mg/dL:** unit for measuring blood glucose
- ⑨ **Battery symbol:** indicates meter battery is running low and needs to be replaced
- ⑩ **Blood insertion symbol:** indicates meter is ready for the application of a drop of blood or control solution
- ⑪ **Pre-meal test flag:** used for tests done before eating
- ⑫ **Post-meal test flag:** used for tests done after eating
- ⑬ **Fasting test flag:** used for tests done after fasting for at least 8 hours
- ⑭ **Month/Day/Hour/Minute**

NOTE

- It is recommended to check if the display screen on the meter matches the illustration above every time the meter turns on.
- **Do not** use the meter if the display screen does not exactly match the illustration as the meter may show incorrect results.

Understanding the Color Screen

Consult your healthcare professional to activate the Target Range Indicator and change the target range for your needs.

The three-color (Red/Green/Blue) screen will display according to the preset target range.

Please see the following diagram to help you understand each color screen.

Target Range Indicator must be turned 'On' to view the blood glucose results on the colored screen. Please see page 27 to activate this function.

- A red screen indicates that the measured value is below the target range, which does not mean hypoglycemia.
- A blue screen indicates that the measured value is above the target range, which does not mean hyperglycemia.

Setting Up the ReliOn™ Premier BLU System

Meter settings, such as time and date, should be checked and updated before using the meter or after changing the meter batteries.


Press and hold the  button for three seconds to enter 'SET' mode. After all settings are finished, press and hold the  button for three seconds to turn off the meter.

Press the < or > button to change values. Press and hold the < or > button to scroll faster.

Step 1 Entering the SET Mode

Press and hold the  button for three seconds. After all the segments flash across the screen, 'SET' is displayed on the screen.

Press the < or > button to select 'YES' and press the  button to go to the next step.

NOTE

- Follow steps 2–3 to pair your meter and smartphone.
- Pairing allows the meter to communicate wirelessly with your smartphone.
- Ensure that devices are within the maximum Bluetooth® range of 32 feet (10 meters).
- Before pairing your meter with your smart device (such as a smartphone), download the compatible app to your smart device.

Setting Up Bluetooth®

Step 2 Bluetooth® Pairing

- ① Download Mobile App:** Download the data management mobile app by visiting relionbgm.com/app and following the instructions to download the compatible mobile app.

- ② Initiate your meter's Bluetooth® pairing mode by one of the following two methods:**

1. On your smartphone, open the app you have downloaded from the section above, titled Download Mobile App.
2. With your meter, take a blood glucose reading. Upon removal of the test strip, the meter will automatically go into Bluetooth® pairing mode.

Or, press and hold down your meter's right arrow button for three seconds. This will place your meter into Bluetooth® pairing mode.

The Bluetooth® (Bluetooth) will blink, indicating your meter is in Bluetooth® pairing mode.

3. Then follow the directions on the smartphone's app screen to complete the pairing.
4. If the smartphone displays 'Unable to Sync', follow the steps listed on the app screen.

NOTE

- The  symbol will appear on the screen when the Bluetooth® feature is on.
- When the  button when 'On' or 'OFF' blinks on the screen.
- Some smartphones, especially those that are not tested or approved by the manufacturer may **NOT** be compatible with your meter.

Adjusting the Date and Time

Step 3 Setting the Year

Press the < or > button until the correct year appears. After setting the year, press the  button to confirm your selection and to go to the next step.

Step 4 Setting the Month

A number indicating the month will blink on the screen.

Press the < or > button until the correct month appears. Press the  button to confirm your selection and to go to the next step.

Step 5 Setting the Date

Press the < or > button until the blinking number shows the correct date. Press the  button to confirm the date and to go to the next step.

Step 6 Setting the Time Format

The meter can be set in the 12-hour (AM/PM) or the 24-hour clock format.

Press the < or > button to select a format. The AM•PM symbol is not displayed in the 24-hour format. After selecting the format, press the $\textcircled{1}$ button to go to the next step.

Step 7 Setting the Hour

Press the < or > button until the correct hour appears. Press the $\textcircled{1}$ button to confirm the hour and go to the next step.

Step 8 Setting the Minute

Press the < or > button until the correct minute appears. After setting the minute, press the $\textcircled{1}$ button to go to the next step.

Turning the Sound On/OFF

Step 9

On pressing the < or > button, the screen will display 'On' or 'OFF'. Press the ① button to confirm the selection.

The meter will beep in the following instances if set to 'On'.

- When a button is pressed to turn on the meter
- When the test strip is inserted in the meter
- When the blood sample is absorbed by the test strip and the test starts
- When the test result is displayed
- When the < button is pressed to set the post-meal (PP2) alarm
- When it is time for a preset blood glucose test

If the sound is set to 'OFF', none of the sound functions will work.

After setting the sound, press the ① button to go to the next step.

NOTE

- The speaker symbol is displayed only when the sound is set to 'OFF'.

Turning on the Test Strip Expiration Date Indicator

Step 10

This setting allows you to turn the test strip expiration date indicator 'On' or 'OFF'. This setting turns the function 'On' or 'OFF' **only**. See page 30 to set the test strip expiration date.

When 'EP' appears on the screen, press the < or > button. The screen will display 'On' or 'OFF'. Press the $\textcircled{1}$ button to confirm the setting and go to the next step.

If you do not want to set the indicator, press the $\textcircled{1}$ button to go to the next step while the screen displays 'OFF'.

NOTE

- If the preset expiration date expires, the meter will display 'EP' when the test strip is inserted.
- The display will alternate between 'EP', and date and time when the test result is displayed right after the test.
- If the expiration date is set to October of 2025, the meter will display EP at the start of November, 2025.

Setting the Target Range

Step 11 Activating the Target Range

Consult your healthcare professional to activate the Target Range Indicator function and change the target range for your needs. Press the < or > button to navigate to the Target Range setting.

Press the  to confirm activating the target range.

Target Range 'OFF'

Target Range 'ON'

Step 12 Setting the Low-end of the Target Range

Press the < or > button to adjust the lowend value (tr Lo) and confirm by pressing the .

Step 13 Setting the High-end of the Target Range

Press the < or > button to adjust the highend value (tr Hi) and confirm by pressing the .

Press the  for three seconds to save the setting and turn the meter off.

NOTE

- You can change the low-end target and the high-end target within the given ranges below:
 - tr Lo: 60-110 mg/dL
 - tr Hi: 90-300 mg/dL
- The low-end of the target range can be changed first. The high-end can then be set to a number that ensures that the user has at least 20 mg/dL within the target range.
 - For example, if the [Tr Lo] is set to 99 mg/dL, the [Tr Hi] can be set to a number in between 121 to 300 mg/dL.
- The color screen will display according to the target range as follows: (the image below is an example only. The color screens will be set according to your target range setting values)

[Below the target range]

[Within the target range]

[Above the target range]

- Color screen will not be displayed when testing with the control solution flag on. Please refer to page 32 step 2 for control solution testing.
- The values of target range indicator can only be set on the meter and are not synchronized with the mobile application.

IMPORTANT

- The color screen helps you to see if your blood glucose results are within your target range. Please note that it should not be used as the sole basis for making clinical decisions.
- A red screen indicates that the measured value is below the target range, which does not mean hypoglycemia.
- A blue screen indicates that the measured value is above the target range, which does not mean hyperglycemia.

Target Blood Glucose Ranges

Reminders

Time of day

*Your target ranges
from your healthcare professional*

Before breakfast

Before lunch or dinner

One hour after meals

Two hours after meals

Expected Values: Normal blood glucose levels for an adult without diabetes are below 100 mg/dL before meals and fasting* and are less than 140 mg/dL two hours after meals.¹

*Fasting is defined as no caloric intake for at least eight hours.

Reference

¹ American Diabetes Association (Standards of Medical Care in Diabetes - 2021. Diabetes Care), January 2021, vol 44 (Supplement 1): S15-S33.

Setting the Test Strip Expiration Date Indicator

Step 1 Entering the Expiration Date Setting

Press and hold the < and > buttons at the same time for three seconds to enter the expiration date settings. After all segments flash across the screen, 'EP' will appear on the display.

NOTE

- The test strip expiration date is printed on the test strip vial.

Step 2 Setting the Expiration Year

A number indicating the year will blink in the left corner of the screen. Press the < or > button until the correct year appears.

Press the  button to confirm the year and go to the next step.

Step 3 Setting the Expiration Month

A number indicating the month will blink at the bottom of the screen. Press the < or > button until the correct month appears.

After setting the month, press the  button for three seconds to save your setting and turn off the meter.

Checking the System

ReliOn™ Premier BLU Blood Glucose Meter and Test Strips should be checked using ReliOn™ Premier Control Solutions, available in two levels (Level 1 and Level 2). The Control Solutions contain a known amount of glucose and are used to check that the meter and test strips are working properly.

The test strip vials have the Control Solution ranges printed on their labels. Compare the result displayed on the meter to the control solution range printed on the test strip vial.

Before using a new meter or a new vial of test strips, you should conduct a control solution test following the procedure on pages 32-34 with ReliOn™ Premier Control Solutions.

NOTE

- Use ReliOn™ Premier Control Solutions only (available for purchase separately).
- Check the expiration date printed on the bottle. When the control solution bottle is first opened record the discard date (date opened plus three months) in the space provided on the label.
- Make sure the meter, test strips, and control solutions are at room temperature before testing. Control solution tests must be done at room temperature, 68–77 °F (20–25 °C).
- Before using the control solution, shake the bottle, discard the first one or two drops, and wipe the tip clean.
- Close the control solution bottle tightly and store at a temperature between 46–86 °F (8–30 °C).

Use ReliOn™ Premier Control Solution when:

- You want to practice the test procedure using the control solution, instead of blood
- You use the meter for the first time
- You begin using a new vial of test strips
- You suspect that the meter or test strips are not working properly
- You think your test results are inaccurate or they do not reflect how you feel
- The meter is dropped or damaged

Control Solution Testing

Step 1

Insert a test strip into the meter's test strip port with the contact bars facing upwards.

Gently push the test strip into the port until the meter beeps. Be careful not to bend the test strip while pushing it in.

The  symbol will show up.

Step 2


You can flag the control solution test result by pressing the > button for three seconds. To undo the control solution flag, press the > button for three seconds again.

Step 3

Shake the ReliOn™ Premier Control Solution bottle before each test.

Remove the cap and squeeze the bottle to discard the first one or two drops.

Wipe the tip with a clean tissue or cloth. Dispense a drop of control solution onto a clean non-absorbent surface. It helps to squeeze a drop onto the top of the cap as shown.

After the symbol appears on the display, apply the solution to the narrow edge of the test strip until the meter beeps. Make sure the confirmation window fills completely.


NOTE

- The meter may switch off if the control solution sample is not applied within two minutes of the symbol appearing on the screen.
- If the meter turns off, remove the test strip, reinsert, and start from step 2.

Step 4

The display segments will rotate clockwise and a test result will appear after the meter counts down from '5' to '1'.

When flagged, the result is stored in the meter's memory but it is not included in the averages.

Step 5

Compare the result displayed on the meter to the range printed on the test strip vial. The result should fall within the range.

CAUTION

- The range printed on the test strip vial is for the ReliOn™ Premier Control Solutions **only**. It has nothing to do with your blood glucose level.

NOTE

- ReliOn™ Premier Control Solutions can be ordered separately by contacting Customer Service at 1-833-256-3561.

Comparing the Control Solution Test Results


The test result of each control solution should be within the range printed on the label of the test strip vial. Repeat the control solution test if the test result falls outside of the range. Out-of-range results may occur in the following situations:

Situations	Do This
<ul style="list-style-type: none">When the control solution bottle was not shaken well.When the meter, test strip or control solution was exposed to high or low temperatures.When the first one or two drops of control solution was not discarded or the tip of the bottle was not wiped clean.When the meter is not functioning properly.	Repeat the control solution test by referring to the NOTE on page 31.
<ul style="list-style-type: none">When the control solution is past the expiration date printed on the bottle.When the control solution is past its discard date (the date the bottle was opened plus three months).When the control solution is contaminated.	Discard the used control solution and repeat the test using a new bottle of control solution.

If results continue to fall outside the range printed on the test strip vial, the ReliOn™ Premier Test Strips and ReliOn™ Premier BLU Blood Glucose Meter may not be working properly. **Do not** use your system and contact Customer Service at 1-833-256-3561.

Using the Lancing Device (sold separately)

A lancing device, such as the ReliOn™ Lancing Device (sold separately), is needed in order to collect a blood sample.

The lancing device is for use by a single user only and should not be shared with anyone including other family members.

Use a soft cloth or tissue to wipe the lancing device. If necessary, a small amount of alcohol on a soft cloth or tissue may be used.

CAUTION

To reduce the chance of infection

- All components that come into contact with blood samples should be considered to be biohazardous and capable of transmitting viral diseases, even after disinfection.
- Do not use a lancet that has been used by another person and never share your meter and lancing device with anyone.
- Wash your hands thoroughly with soap and water after handling your meter, lancing device or test strips.
- Keep your meter and lancing device clean. See the cleaning and disinfection section on page 54.

NOTE

- Repeated puncturing at the same sample site may cause pain or skin calluses (thick hard skin).
- Choose a different site each time you test.
- The ReliOn™ Lancing Device is given only as an example. Users should carefully read and follow the instructions for use provided with the lancing device they have purchased.

Preparing the Lancing Device

Step 1


Wash hands and sample site with soap and warm water.

Rinse and dry thoroughly.


Step 2

Pull the adjustable cap sideways to remove.


Step 3

Insert a new disposable lancet firmly into the lancet holder. Twist off the protective cover of the lancet and set it aside, then replace the adjustable tip. Keep the protective cover to replace on top of the used lancet after testing.


Step 4

The lancing device has five puncture depth settings. The shorter lines are for a shallower puncture and the longer lines are for a deeper puncture. Adjust the depth of puncture as necessary.

Step 5

Cock the lancing device by holding the body in one hand and pulling on the charger. A click may be heard.

NOTE


- The skin depth to retrieve samples will vary for various people at different sample sites.
- The lancing device's adjustable tip allows the best depth of skin penetration for an adequate sample size.
- A beginning setting of '3' is recommended.

Preparing the Meter and Test Strip

Step 6

Insert a test strip with the contact bars facing upwards into the meter's test strip port.


Push the test strip in gently until the meter beeps. Be careful not to bend the test strip. The symbol will appear on the screen.


Applying Blood Sample

Step 7

Obtain a blood sample using the lancing device. Place the device against the pad of the finger. The best puncture sites are on the middle or ring fingers.

Press the release button. Remove the device from the finger. Wait a few seconds for a blood drop to form.


A minimum volume of 0.5 µL is needed to fill the confirmation window.

Step 8

After the symbol appears on the screen, apply the blood sample to the narrow end of the test strip until the meter beeps.

If the confirmation window is not filled in time because of abnormal viscosity (thickness and stickiness) or insufficient volume, the **Er4** message may appear.

It is recommended to place the test strip vertically into the blood sample site as shown below.

NOTE

- The meter may switch off if the blood sample is not applied within two minutes of the  symbol appearing on the screen.
- If the meter turns off, remove the test strip, reinsert it and start from step 2.

Step 9

At this time, the display segments will rotate clockwise while the blood sample is going into the test strip.


The test result will appear after the meter counts down from '5' to '1'.

The result will be automatically stored in the meter's memory. If the test result is within the preset target value, the result appears with the green-colored backlight.

If the test strip is removed after the test result is displayed, the meter will automatically switch off after three seconds

Discard used test strips safely in disposable containers.

If the Bluetooth® feature is activated, the meter will send the test result to the connected smart device (such as smartphone).

NOTE

- If the test result is lower than the target value, the meter displays result with the red screen.
- If the test result is higher than the target value, the meter displays result with the blue screen.

NOTE

To transmit glucose data using the Bluetooth® feature:

- The Bluetooth® feature on the meter must be turned on
 - The meter and a smart device (such as a smartphone) must be paired
 - The paired smartphone app must be launched
- The meter will transmit data in the following cases:
- When the test strip is ejected after measuring
 - When the meter is turned on (when untransmitted data exists)

Step 10

You can attach a flag to a result to indicate particular situations while the test strip is still in the meter.

When the result is displayed right after a test, press the < or > button to select a pre-meal flag (🍎), a post-meal flag (🏃), a fasting flag (🚶) or a control solution flag (⌚).

When you remove the test strip while the desired flag is blinking, the test result is stored with the flag.

If you do not want to add any flags on the test result, remove the test strip after the test result is displayed.

No flag

Pre-meal flag

Post-meal flag

Fasting flag

Control solution flag


NOTE

The colored backlight display will not appear when attaching a flag on the test result.

Discarding Used Lancets

Step 1


Pull the adjustable cap sideways to remove.

Step 2

Place the protective cover on the lancet.

Remove lancet and dispose of the used lancet in a proper biohazard container.

CAUTION


- Check for damage before using the lancet. If it has been damaged, please discard it and use another lancet.
- The lancet is very sharp. Please keep away from children.
- Keep the lancets in a cool and dry place.
- The lancet is for single use only. Never share or reuse a lancet. Always dispose of lancets properly.

'HI' and 'Lo' Messages

'HI' Message

The meter displays results between 20–600 mg/dL.

'HI' appears when the blood glucose level is greater than 600 mg/dL and indicates severe hyperglycemia (much higher than normal glucose levels).

If 'HI' is displayed again upon retesting, please contact your healthcare professional immediately.

'Lo' Message

'Lo' appears when a test result is less than 20 mg/dL and indicates severe hypoglycemia (very low glucose levels).

If 'Lo' is displayed again upon retesting, please contact your healthcare professional immediately.

NOTE

- When the meter displays 'HI' or 'Lo', the color screen will not be displayed.
- If the messages persist, please contact Customer Service at 1-833-256-3561. At all other times or in case of emergency, contact your healthcare professional or emergency medical response.

Transferring Test Results Using the Data Cable

Test results stored in ReliOn™ Premier BLU Blood Glucose Meter can be transferred from the meter to a compatible desktop data management application via cable.

'PC' is displayed when the data cable connects the meter with a computer.

For more information, visit relionbgm.com/app or contact Customer Service at 1-833-256-3561.

NOTE

- Test results can also be transferred to paired smart device app wirelessly using Bluetooth®.
- To pair your meter and smart device, see page 21.

Meter Memory

The ReliOn™ Premier BLU Blood Glucose Meter can save up to 1,000 blood glucose test results with time and date. If the memory is full, the oldest test result will be deleted and the latest test result will be stored.

The meter calculates and displays the averages of total test results, pre-meal (🍎) test results, post-meal test (🥗) and fasting test results (🍽️) from the last 1, 7, 14, 30 and 90 days.

Viewing Past Test Averages

Step 1

Press the <, > or ⓘ button to turn on the meter.


The current date and time will be displayed at the bottom of the screen followed by the one-day average value and the number of the test results saved within the current day.

The number of tests within the current day


Step 2 Viewing Averages

Press the < button to view 7-, 14-, 30- and 90-day average values, and the number of tests performed for the last test period.


Step 3 Viewing Pre-meal Averages

Continue to press the < button to view 1-, 7-, 14-, 30- and 90-day average values and the number of tests performed pre-meals with the  symbol for the last test period.


Step 4 Viewing Post-meal Averages

Press the < button to view 1-, 7-, 14-, 30- and 90-day average values and the number of tests performed post-meals with the  symbol for the last test period.

Step 5 Viewing Fasting Averages

Press the < button to view 1-, 7-, 14-, 30- and 90-day average values and the number of tests performed during fasting with the  symbol for the last test period.

Step 6 Viewing Pre-meal Averages

Use the > button to scroll through the averages seen previously.

Press the ⓧ button to turn off the meter.

NOTE

- The control solution test results saved with the  symbol are **not included** in the averages.

Viewing Test Results

Step 1

Press the <, > or ① button to turn the meter on.


The current date and time will be displayed on the bottom of the screen followed by the one-day average value and the number of test results saved within the current day.

The number of tests within the current day

Step 2

Use the > button to scroll through the test results, starting from the most recent and ending with the oldest.

Press the < button to return to the result seen previously.

After checking the stored test results, press and hold the ① button to turn off the meter.

NOTE

- The control solution test results saved with the  symbol will be displayed with  symbol when you review the stored test results.

Setting the Alarm Function

Four types of alarms can be set in the ReliOn™ Premier BLU Blood Glucose Meter. One post-meal alarm (PP2 alarm) and three time set alarms (alarm 1–3).


- The PP2 alarm goes off two hours after setting the alarm.
- The alarms ring for 15 seconds and can be silenced by pressing the <, > or \odot button, or by inserting a test strip.

Setting the Post-meal Alarm (PP2 alarm)

Step 1 Turning the PP2 Alarm On

Without inserting a test strip, press and hold the < button for three seconds to set the post-meal alarm.

'PP2', the  symbol and 'On' will be displayed. The screen will then automatically change to the memory recall mode.

At this time, the  symbol, indicating that the PP2 alarm has been set, will be displayed on the screen.


NOTE

- The PP2 alarm will automatically turn off if the meter's time setting is adjusted to more than two hours before or just past the currently activated PP2 alarm time.

Step 2 Turning the PP2 Alarm OFF

To turn off the PP2 alarm, press and hold the < button for three seconds. 'PP2', the  symbol and 'OFF' will appear on the screen.

Then the screen will change automatically to the memory recall mode without the  symbol displayed.

Setting the Time Alarms (Alarm 1–3)

Step 1

Without inserting a test strip, press the < and Ⓛ buttons simultaneously for three seconds to enter the time alarm setting.

'alarm 1' will be displayed while 'OFF' blinks on the screen.

Step 2

On pressing the > button, 'alarm 1' is set and 'On' is displayed on the screen. Press the > button again to cancel 'alarm 1'. 'OFF' will blink on the screen.

Step 3

Press the < button to adjust the time of 'alarm 1'. A number indicating the hour will blink on the screen. Press the > button to set the hour.

Step 4

On pressing the < button, the number indicating the minute will start blinking.
Press the > button to set the minute.

Step 5

Press the ① button to finish and to go to 'alarm 2' setting.

Repeat steps 2 to 5 to set the remaining time alarms (alarm 2–3).

Step 6

Press the ① button for three seconds to finish and turn the meter off.

Caring for the ReliOn™ Premier BLU System

To minimize the risk of transmission of blood-borne pathogens, the cleaning and disinfection procedure should be performed as recommended in the instructions below.

Wash your hands thoroughly with soap and water after handling the meter, lancing device, or test strips.

If the meter is being operated by a second person who is providing testing assistance to the user, the meter and lancing device should be disinfected prior to use by the second person.

Cleaning and Disinfection

The cleaning procedure is needed to clean dirt as well as blood and other body fluids on the exterior of the meter and lancing device before performing the disinfection procedure.

The disinfection procedure is needed to prevent transmission of blood-borne pathogens.

For the meter and lancing device, this cleaning and disinfection procedure should be performed **once per week**.

NOTE

- The life span of a ReliOn™ Premier BLU Blood Glucose Meter is five years. We recommend disinfecting both the meter and lancing device at least once per week. We have validated a total of 260 cleaning and disinfecting cycles (260 cleaning and 260 disinfection cycles) to represent weekly cleaning and disinfecting over the use life of your meter and lancing device.
- 1 cleaning and 1 disinfection cycle per week
 $X 52 \text{ weeks per year} X 5 \text{ years} = 260 \text{ cleaning and 260 disinfection cycles.}$

We have validated **Clorox Healthcare Bleach Germicidal Wipes** with 0.55 % sodium hypochlorite as the active ingredient for disinfecting the ReliOn™ Premier BLU Meter and lancing device. It has been shown to be safe for use with the meter and lancing device.

This disinfectant is available commercially in towelette form. In addition to ReliOn™ Premier BLU Blood Glucose Monitoring System instructions, **please read the instructions provided by the manufacturer of Clorox Healthcare Bleach Germicidal Wipes before using it.**

Name	Clorox Healthcare® Bleach Germicidal Wipes	
Manufacturer	Clorox® Professional Products Company Phone: 1 800 492 9729 Website: www.cloroxpro.com	
EPA registration number	67619-12	
Active ingredients	Sodium Hypochlorite: 0.55 %	

NOTE

- The disinfectant products can be purchased through online retailers (e.g. Walmart.com) or by telephone at the manufacturer.
- To find out where to purchase the disinfectant product, please contact the manufacturer or visit their website as listed above.

Cleaning and Disinfection Procedures

- ① Open the cap of the Clorox Healthcare Bleach Germicidal Wipes container and pull out one towelette and close the cap.

- ② Wipe the entire surface of the meter three times horizontally and three times vertically using one towelette to pre-clean blood and other body fluids.

- ③ Dispose of the used towelette in a trash bin.

- ④ Pull out one new towelette and wipe the entire surface of the meter three times horizontally and three times vertically using a new towelette to disinfect the meter.

- ⑤ Dispose of the used towelette in a trash bin.

- ⑥ Allow exteriors to remain wet for one minute.

- ⑦ Repeat the same procedure for the lancing device (step ① to step ⑥).

NOTE

- After the cleaning and disinfection procedure, the control solution should be tested to confirm that the meter works properly before using the meter.
- Control solution tests should be performed using two different levels of ReliOn™ Premier Control Solutions (Level 1 and/or Level 2).
- Verify that the test results are within the range printed on the test strip vial.
- See page 32-34 for how to perform a control solution test.

NOTE

If any of the following deterioration signs appear after cleaning or disinfecting, please stop using the system and contact Customer Service at 1-833-256-3561.

- The inscriptions on the exterior of the meter or lancing device have been removed.
- The color of the meter or lancing device has changed or faded.
- Cracks or roughness develop on the meter or lancing device.
- A part of the segment on the meter display does not flash.
- Control solution test results do not fall within the stated range on the test strip vial.

CAUTION

- **Do not** use other cleaners or disinfectants because other chemicals have not been validated and may damage the meter.
- **Do not** get fluids inside the meter through the test strip port, data transmission port or battery compartment. Never immerse the meter or hold it under running water because this will damage the meter.

Storage and Handling Caution

- Store the meter in a cool and dry place between 32-122 °F.
- **Do not** expose the meter to direct sunlight, heat or excessive humidity for an extended period of time. It is recommended to store and use the test system indoors.
- **Do not** let dirt, dust, blood, or water enter into the meter's test strip port.
- **Do not** drop the meter or subject it to strong shock.
- **Do not** try to fix or alter the meter in any way.
- Store all meter components in the carrying case to prevent loss and help keep the meter clean.
- Avoid getting any liquid or moisture in the test strip vial. This can affect the test strips and cause inaccurate test results.
- **Do not** apply samples other than capillary whole blood or control solution to the test strip.

NOTE

- For additional information or technical assistance contact Customer Service at 1-833-256-3561.

CAUTION

- **Do not** allow any foreign substances or liquid substances, such as dirt, blood, or water, enter into the meter. The meter may be damaged or may malfunction. Follow the warning information provided below to prevent possible damage to the meter.
- **Do not** apply the blood or control solution samples directly to the test strip port.
- **Do not** apply the blood or control solution samples to the test strip while holding the meter in a way that the tip of the test strip faces upwards. The blood or control solution samples may run down the surface of the test strip and flow into the test strip port.
- **Do not** store your meter in unsanitary or contaminated sites.
- Make sure to follow the Cleaning and Disinfection Procedures found in the Caring for the ReliOn™ Premier BLU System section of this user manual.

Understanding Error Messages

Message	What it Means	What to Do
	A used test strip was inserted.	Repeat the test with a new test strip.
	The blood or control solution sample was applied before the  symbol appeared.	Repeat the test with a new test strip and wait until the  symbol appears before applying the blood or control solution sample.
	The temperature during the test was above or below the operating range.	Move to an area where the temperature is within the operating range 42.8–111.2 °F (6–44 °C). Repeat the test after the meter and test strips have reached a temperature within the operating range.
	The blood sample is of insufficient volume or blood is too thick.	Repeat the test after inserting a new test strip.

Message	What it Means	What to Do
	This error message may appear when the wrong blood glucose test strip is used instead of a ReliOn™ Premier Blood Glucose Test Strip.	Repeat the test with a ReliOn™ Premier Blood Glucose Test Strip.
	There is a problem with the meter.	Do not use the meter. Contact Customer Service: 1-833-256-3561.
	There is a problem with the Bluetooth® communication.	Contact Customer Service: 1-833-256-3561.
	An electronic error occurred during the test.	Repeat the test with a new test strip. If the error message persists, contact Customer Service: 1-833-256-3561.

NOTE

For questions or concerns, contact Customer Service at 1-833-256-3561 seven days a week. Standard hours are Monday-Friday 7a.m. to 7 p.m. Central time. Saturday-Sunday (and holiday) hours are 8 a.m. to 5 p.m. Central time.

General Troubleshooting

Problem	Troubleshooting
The display is blank even after inserting a test strip.	<ul style="list-style-type: none">Check whether the test strip is inserted with the contact bars facing up. Check if the test strip has been inserted completely into the test strip port.Check if the appropriate test strip was used.Check whether the batteries are inserted with the '+' side facing up.Replace the batteries.
The test does not start even after applying the blood sample on the test strip.	<ul style="list-style-type: none">Check if the confirmation window is filled completely.Repeat the test after inserting a new test strip.
The test result does not match the way you feel.	<ul style="list-style-type: none">Repeat the test after inserting a new test strip.Check the expiration date of the test strip.Perform control solution test.

NOTE

For questions or concerns, contact Customer Service at 1-833-256-3561 seven days a week. Standard hours are Monday-Friday 7a.m. to 7 p.m. Central time. Saturday-Sunday (and holiday) hours are 8 a.m. to 5 p.m. Central time. At all other times or in case of emergency, please contact your healthcare professional or emergency medical response.

Performance Characteristics

Performance of the ReliOn™ Premier BLU Blood Glucose Monitoring System has been evaluated in laboratory and clinical tests.

Accuracy

The accuracy of the ReliOn™ Premier BLU Blood Glucose Monitoring System was assessed by comparing blood glucose results obtained by patients with those obtained using a YSI Model 2300 Glucose Analyzer, a laboratory instrument. The following results were obtained by diabetic patients at clinic centers.

Slope	1.0223
Y-Intercept	-1.3686
Correlation coefficient (r)	0.9934
Number of subjects	371
Range tested	48–553 mg/dL

Accuracy results for glucose concentration < 75 mg/dL

Within ± 5 mg/dL	Within ± 10 mg/dL	Within ± 15 mg/dL
25/41 (61 %)	40/41 (97.6 %)	41/41 (100 %)

Accuracy results for glucose concentration ≥ 75 mg/dL

Within ± 5 %	Within ± 10 %	Within ± 15 %	Within ± 20 %
231/330 (70 %)	317/330 (96.1 %)	330/330 (100 %)	330/330 (100 %)

User performance results for glucose concentrations between 48 mg/dL and 553 mg/dL.

Within \pm 15 mg/dL and Within \pm 20 %
371/371 (100 %)

Precision

Precision studies were performed in a laboratory using the ReliOn™ Premier BLU Blood Glucose Monitoring System.

Within Run Precision			
Blood average	43 mg/dL	SD = 1.9 mg/dL	CV = 4.3 %
	71 mg/dL	SD = 2.1 mg/dL	CV = 2.9 %
	135 mg/dL	SD = 3.8 mg/dL	CV = 2.8 %
	203 mg/dL	SD = 5.2 mg/dL	CV = 2.6 %
	343 mg/dL	SD = 11 mg/dL	CV = 3.2 %

Between Run Precision			
Control average	36 mg/dL	SD = 1.4 mg/dL	CV = 3.8 %
	114 mg/dL	SD = 3.4 mg/dL	CV = 3.0 %
	135 mg/dL	SD = 8.2 mg/dL	CV = 2.4 %

EMC Table

The following tables contain the Manufacturer's declaration and additional information required by IEC 60601-1-2:2014+AMD1:2020

Guidance and Manufacturer's Declaration – Electromagnetic Emission		
Emission test	Compliance	Electromagnetic environment - guidance
RF Emission CISPR 11	Group 1	The meter uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emission CISPR 11	Class B	
Harmonic Emission IEC 61000-3-2	N/A	
Voltage fluctuation and flicker IEC 61000-3-3	N/A	The meter is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Guidance and Manufacturer's Declaration – Electromagnetic Emission

The ReliOn™ Premier BLU is intended for use in the electromagnetic environment specified below.

The customer or the user of the ReliOn™ Premier BLU system should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, a relative humidity of at least 30% is recommended.
Electrical fast transient/burst IEC 61000-4-4	±2 kV (for power supply lines) ±1 kV (for input/output lines)	±2 kV (for power supply lines) ±1 kV (for input/output lines)	Mains power quality should be that of a typical commercial or hospital environment.
Surge IED 61000-4-5	±0.5 kV, ±1 kV differential mode ±0.5 kV, ±1 kV, ±2 kV common mode	±0.5 kV, ±1 kV differential mode ±0.5 kV, ±1 kV, ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.

Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical commercial or hospital environment.
Proximity magnetic fields in the frequency range 9 kHz to 13.56MHz immunity IEC 61000-4-39	8 A/m 30kHz CW Modulation 65 A/m 134.2 kHz PM 2.1 kHz 7.5 A/m 13.56 MHz PM 50kHz	8 A/m 30kHz CW Modulation 65 A/m 134.2 kHz PM 2.1 kHz 7.5 A/m 13.56 MHz PM 50kHz	Resistance to magnetic fields was tested and applied only to surfaces of enclosures or accessories accessible during intended use
NOTE: U_T is the main voltage (AC) prior to the application of the test level.			
a For example, 10/12 means 10 cycles at 50 Hz or 12 cycles at 60 Hz.			

Recommended separation distances between portable and mobile communication equipment and the meter system

The ReliOn™ Premier BLU is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the meter as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance according to frequency of transmitter [m]	
	IEC 60601-1-2: 2014	
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 2.7 GHz $d = 2.0\sqrt{P}$
0.01	0.12	0.20
0.1	0.38	0.63
1	1.2	2.0
10	3.8	6.3
100	12	20

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1) At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2) These guidelines may not apply in all situations.

Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Guidance and manufacturer's declaration - electromagnetic immunity			
The meter is intended for use in the electromagnetic environment specified below. The customer or the user should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment guidance
Conducted RF IEC 61000-4-6	<p>3 Vrms 150 kHz – 80 MHz</p> <p>Outside ISM Bands amateur radio bands Bandsd</p> <p>6 Vrms 150 kHz – 80 MHz In ISM bands amateur radio bands Bandsd</p>	<p>3 Vrms</p> <p>6 V/m</p>	<p>Portable and mobile RF communications equipment should not be used closer to any part of the Ultrasound System, including cables, than the recommended separation distance. This is calculated using the equation applicable to the frequency of the transmitter.</p> <p>Recommended Separation Distance $d = 1.2/\sqrt{P}$</p> <p>IEC 60601-1-2:2014 $d=2.0 \text{ m}$ 80 MHz to 2.7 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of marked with the following</p> 
Radiated RF IEC61000-4-3	3 V/m 80 MHz to 2.7 GHz		

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations, Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcasts and TV broadcasts cannot be predicted theoretically with accuracy.

To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the meter is used exceeds the applicable RF compliance level above, the meter should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the ReliOn™ Premier BLU meter.

b In the frequency range exceeding 150 kHz – 80 MHz, the electric field strength should be not higher than 3 V/m.

c The ISM (Industrial, Scientific and Medical) bands between 150 kHz and 80 MHz are 6.765 MHz to 6.795MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz.

d The amateur radio bands between 0.15 MHz and 80 MHz are 1.8 MHz to 2.0 MHz, 3.5 MHz to 4.0 MHz, 5.3 MHz to 5.4 MHz, 7 MHz to 7.3 MHz, 10.1 MHz to 10.15 MHz, 14 MHz to 14.2 MHz, 18.07 MHz to 18.17 MHz.

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

This product is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. Portable RF communications equipment should be used no closer than 30 cm (12 inches) to any parts of the product. Otherwise, the performance of this equipment could be degraded.

Guidance and manufacturer's declaration - electromagnetic immunity					
The meter is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. Portable RF communications equipment should be used no closer than 30cm (12 inches) to any part of the meter. Otherwise, degradation of the performance of this equipment could result.					
Immunity test	Band ^a	Service ^a	Modulation	IEC60601 test level	Compliance level
Proximity fields from RF wireless Communications IEC61000-4-3	380–390 MHz	TETRA 400	Pulse modulation 18Hz	27 V/m	27 V/m
	430–470 MHz	GMRS 460 FRS 460	FM ±5 kHz deviation 1 kHz sine	28 V/m	28 V/m
	704–787 MHz	LTE Band13, 17	Pulse modulation 217 Hz	9 V/m	9 V/m
	800–960 MHz	GSM800:900 TETRA 800 iDEN 820 CDMA 850 LTE Band 5	Pulse modulation 18Hz	28 V/m	28 V/m
	1700–1990 MHz	GSM 1800 CDMA 1900 GSM 1900 DECT LTE Band 1,3,4,25 UMTS	Pulse modulation 217 Hz	28 V/m	28 V/m
	2400–2570 MHz	Bluetooth WLAN 802.11b/g/n RFID 2450 LTE Band 7	Pulse modulation 217 Hz	28 V/m	28 V/m
	5100–5800 MHz	WLAN 802.11a/n	Pulse modulation 217 Hz	9 V/m	9 V/m

NOTE: If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1m. The 1m test distance is permitted by IEC 61000-4-3.

- a** For some services, only the uplink frequencies are included.
- b** The carrier shall be modulated using a 50% duty cycle square wave signal.
- c** As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be the worst case.

Warranty Information

Warranty

The manufacturer warrants that the ReliOn™ Premier BLU Blood Glucose Meter shall be free of defects in material and workmanship in normal use for a period of five years.

The meter must have been subjected to normal use. The warranty does not cover improper handling, tampering, use or service of the meter.

Any claim must be made within the warranty period.

The manufacturer will, at its discretion, repair or replace a defective meter, or meter part that is covered by this warranty.

As a matter of warranty policy, the manufacturer will not reimburse the consumer's purchase price.

Obtaining Warranty Service

To obtain warranty service, you must return the defective meter or meter part along with proof of purchase.

Returns

For instructions on how to return your meter, contact Customer Service at 1-833-256-3561.

Meters returned without this authorization will not be accepted.

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